

LIFE CAN BE BUSY AND REMEMBERING TO TAKE MEDICATION IS NOT ALWAYS EASY. **HEALTHBEACON** IS HERE TO SUPPORT YOU TO STAY ON TRACK WITH YOUR MEDICATION SCHEDULE.



WHAT DOES HEALTHBEACON DO?

REMINDS

HealthBeacon will light up on the day of your scheduled injection **customized** to your preferred time. You will also receive a text message indicating that it is time to take your medication.



TRACKS

The HealthBeacon helps you track your progress by calculating and displaying a **Personal Adherence Score (PAS)**. This score is displayed on your HealthBeacon screen and represents the number of on-time injections you have taken to date.



SAFELY DISPOSES

The HealthBeacon **safely stores** used injections. The inner red sharps bin will need to be replaced when full and a new sharps bin will be automatically sent to you.




SUPPORTS

Our **Customer Care team** is on hand for any questions regarding your HealthBeacon. For example, if you are away from your HealthBeacon and take your medication, just call our Customer Care team at **support.usa@healthbeacon.com** and they will record that you have taken your medication and will update your Personal Adherence Score accordingly.



HEALTHBEACON

INJECTION CARE MANAGEMENT SYSTEM



YOUR RIGHTS AND CHOICES

YOU MAY CONTACT HEALTHBEACON AT ANY TIME TO WITHDRAW YOUR CONSENT FROM THE PROGRAM. YOU MAY CONTACT HEALTHBEACON TO REQUEST ACCESS TO YOUR PERSONAL DATA AND THAT WE CORRECT OR UPDATE ANY OF YOUR PERSONAL DATA. NOTE THAT IF YOU WITHDRAW CONSENT OR REQUEST RESTRICTION OR DELETION OF YOUR PERSONAL DATA, YOU MAY NO LONGER BE ABLE TO PARTICIPATE IN THE PROGRAM OR PARTS OF IT. YOUR PERSONAL DATA IS IMPORTANT TO US. PLEASE CONTACT THE HEALTHBEACON CUSTOMER CARE TEAM AT SUPPORT.USA@HEALTHBEACON.COM FOR MORE INFORMATION.



LET'S GET STARTED!

FOLLOW THESE 3 STEPS AND YOU WILL BE READY TO USE THE HEALTHBEACON.

STEP 1: CONNECT YOUR HEALTHBEACON

First you'll need to place your HealthBeacon in a location with the following:

1. An electric outlet. The HealthBeacon needs to be plugged in at all times.
2. Good cell signal strength.



CONFIRM ACTIVATION

Your pre-programmed personal medication schedule will appear on the screen once you plug in your unit and activation is successful. This may take up to 15 minutes.

SCREEN WILL DISPLAY:

1. The date of your next scheduled dose.
2. Your **Personal Adherence Score** - this will change to 100% after your first on-time injection drop.



Please note: It may take 15 minutes for your HealthBeacon to connect to the network. If your next medication date does not display after 15 minutes please contact our Customer Care team for support.

STEP 2: ON THE DAY OF YOUR INJECTION

1. The **HealthBeacon** will light up showing a blue reminder on top of the display screen.
2. You will receive a text message letting you know that your injection is due.
3. **Dispose** of your injection:
 - Place your used injection in your HealthBeacon and dispose of caps into recycling bin.
 - Wait until you see a flash - this is the HealthBeacon taking a record of your disposed injection.



- Close the door. "Well Done" will appear on the screen.



Your Personal Adherence Score will update shortly and the date of your next injection will be displayed on the screen.



STEP 3: SHARPS BIN REPLACEMENT

BIN ALMOST FULL

The **HealthBeacon** will let you know when the red inner sharps bin is almost full and will automatically send you a new one.



BIN FULL

When the "Bin Full" message appears on the screen its time to **replace your red sharps bin**.



1. **Open** the side door and remove the full sharps bin.
2. **Lock** the sharps bin by sliding the top lid shut.
3. **Slide** the new replacement bin inside the HealthBeacon, ensuring that the door on the white lid is open and on the same side as the front of the HealthBeacon unit.



DO NOT RETURN THE SHARPS BIN IN THE MAIL

Please dispose of the red sharps bin according to your state guidelines. More information regarding state guidelines can be found at: www.SafeNeedleDisposal.org

CUSTOMER CARE TEAM

Contact **HealthBeacon Customer Care** at support.usa@healthbeacon.com if you have any questions regarding:

- Schedule change
- Vacation
- Self reporting a drop
- Bin replacements

TOP TIPS!

KEEP YOUR HEALTHBEACON PLUGGED IN AT ALL TIMES!

PUSH THE SNOOZE BUTTON TO:

1. View your injection site location reminder.
2. To disable the blue lights for one hour.



FREQUENTLY ASKED QUESTIONS

1. IS THE UNIT SAFE TO BE PLUGGED IN ALL THE TIME?

Yes, it is safe and consumes minimal power. The unit must remain plugged in at all times in order to work properly.

2. CAN OTHER MEDICATIONS BE DISPOSED OF IN THE HEALTHBEACON?

No, if you drop other medication in the unit it will be misread as an injection drop and could impact your scheduled reminders and affect your personal adherence score.

3. WHAT SHOULD I DO IF I NEED TO MOVE THE HEALTHBEACON?

First, ensure that the new location has good cell reception. Simply unplug and place your unit in the new location. Remember to plug the unit back in. HealthBeacon will automatically reconnect to the network.

4. WHAT IF I HAVE A QUESTION ABOUT MY TREATMENT?

Please consult with your medical professional directly.

5. WHAT IF THE BIN FULL SCREEN IS STILL APPEARING?

If the Bin Full screen is still appearing after replacing the sharps bin, please notify our Customer Care team.